

# Eliminate Downtime and Calibrate In-House

Calibrate in-house with these easy to use Fluke Calibration models that will help you save money and downtime.

## Temperature Calibration

### Fluke Calibration 9100S/9102S Handheld Dry-Well



Small, light, and portable dry wells.

- Ranges from -10 °C to 375 °C
- Accuracy to  $\pm 0.25$  °C, stability of  $\pm 0.05$  °C at 0 °C

#### Handheld Dry-Well Offer



Receive a FREE Fluke T3000 FC  
Wireless Temperature Module  
valued at **\$371.00**

### Fluke Calibration 9142/9143/9144 Field Metrology Wells



Lightweight and portable dry wells

- Cool to -25 °C in 15 minutes and heat to 660 °C in 15 minutes

**Purchase with a case or insert**

#### Field Metrology Wells Offer



Receive a FREE  
Fluke 5615 or 5609  
Temperature Probe  
valued at **over \$2,000.00**

### Fluke Calibration 9009 Dual-Block Thermometer Calibrator



Double your productivity or cut your calibration time in half.

- Temperatures from -15 °C to 350 °C in one unit
- Two wells in each block for simultaneous comparison calibrations

#### Dual-Block Thermometer Calibrator Offer



Receive a FREE Fluke 179 True RMS  
Multimeter  
valued at **\$653.00**

### Fluke Calibration 6102/7102/7103 Thermometer Calibrators



Portable and extremely stable.

- Calibrates sensors of different size or shape
- Accuracy  $\pm 0.25$  °C

#### Thermometer Calibrators Offer



Receive a FREE Fluke Calibration  
9310/9311/9317 Carrying Case  
valued at **\$1,652.00**

\*Terms and Conditions: Sponsored by Fluke Australia Pty Ltd. Gift with purchase is available for all purchases of a "Qualifying Fluke Purchase" made from an authorised participating Fluke distributor in Australia and New Zealand between 04th May and 21st August 2020 or until stocks last. Redemption gift will be provided by your Fluke distributor. Gift offerings are "up to" the listed dollar value. Prices are in \$AUD and are Fluke Australia's Recommended Trade Price inclusive of GST. Customers are limited to two redemption claims only. Offers cannot be combined with any other discounts, promotions or other specials offers from Fluke Australia or the participating authorised distributor. Fluke reserves the right to substitute a gift of equivalent value or reject the submission if it doesn't meet the redemption requirements. All trademarks are the property of their respective owners. Specifications subject to change without notice. Allow 6 weeks for gift delivery. Additional terms and conditions apply.

# Eliminate on Downtime and Calibrate In-House

## Electrical Calibration

### Fluke Calibration 2638A Hydra Series III Data Acquisition System/DMM



A price performance breakthrough in standalone data acquisition system

- DC accuracy of 0.0024 %
- Best-in-class thermocouple accuracy of 0.5 °C

### Fluke Calibration 1586A Super-DAQ Precision Temperature Scanner and DAQ-STAQ Multiplexer



Best-in-class temperature measurement accuracy and up to 40 isolated input channels

- PRTs:  $\pm 0.005$  °C; Thermocouples:  $\pm 0.29$  °C; Thermistors:  $\pm 0.002$  °C
- Scan speed of up to 10 channels per second

#### Data Acquisition Calibration Scanner Offer



Receive a **FREE**  
Fluke 2638/CASE Transit Case  
**valued at \$1,698.00**

#### Data Acquisition Calibration Scanner Offer



Receive a **FREE**  
Fluke 1586-CASE Carrying Case or  
Fluke 1586/DS-CASE Carrying Case  
**valued at \$1,698.00**

## Pressure Calibration

### Fluke Calibration 700HPPK Pneumatic Test Pump Kit



Rugged and portable way to generate pressure in the field quickly, safely and easily.

- Generates and adjusts pneumatic pressure up to 21 MPa
- Reaches pressure in 20 seconds to full scale into a 30 cm<sup>3</sup> volume

#### Pneumatic Test Pump Kit Offer



Receive a **FREE** Fluke 700G29  
Pressure Gauge Calibrator  
**valued at \$1,360.00**

For the full list of models applicable to this promotion, please reach out to your distributor.

#### How to redeem your **FREE** gift:

1. Purchase a Fluke test tool within this promotion between 04 May - 21 August 2020.
2. Advise your Fluke Authorised Distributor which free Fluke gift you select
3. Your selected free Fluke gift will be provided by your distributor

\*Terms and Conditions: Sponsored by Fluke Australia Pty Ltd. Gift with purchase is available for all purchases of a "Qualifying Fluke Purchase" made from an authorised participating Fluke distributor in Australia and New Zealand between 04th May and 21st August 2020 or until stocks last. Redemption gift will be provided by your Fluke distributor. Gift offerings are "up to" the listed dollar value. Prices are in \$AUD and are Fluke Australia's Recommended Trade Price inclusive of GST. Customers are limited to two redemption claims only. Offers cannot be combined with any other discounts, promotions or other specials offers from Fluke Australia or the participating authorised distributor. Fluke reserves the right to substitute a gift of equivalent value or reject the submission if it doesn't meet the redemption requirements. All trademarks are the property of their respective owners. Specifications subject to change without notice. Allow 6 weeks for gift delivery. Additional terms and conditions apply.

# Eliminate on Downtime and Calibrate In-House

Purchase this Model	Receive this Free		Save
9100S-A-256	Fluke Calibration T3000 FC Temperature Probe		\$371.00
9100S-B-256			
9100S-C-256			
9100S-D-256			
9100S-P-256			
9100S-R-256			
Purchase this Model	and this	Receive this Free	Save
9142-A-256	Fluke Calibration 914X-INSX <b>or</b> Fluke Calibration 9142-CASE	Fluke Calibration 5615 or 5609 Temperature Probe	over \$2000.00
9143-A-256			
9144-A-256			

Purchase this Model	Receive this Free		Save
9009-B-256	Fluke 179 True RMS Mutltimeter		\$653.00
9009-Y-256			

Purchase this Model	Receive this Free		Save
6102-256	Fluke Calibration 9310, 9311 or 9317 Carry Case		up to \$1652.00
7102-256			
7103-256			

Purchase this Model	Receive this Free		Save
700HPPK-NPT	Fluke Calibration 700G29 Pressure Gauge Calibrator		\$1360.00
700HPPK-BSP			
700HPPK-MET			

Purchase this Model	Receive this Free	Save
2638A/20 220-6	Fluke Calibration 2638 or 1586 Carry Case	up to \$1698.00
2638A/20/C 220		
2638A/20/C 240		
2638A/40 220-6		
2638A/40/C 220		
2638A/40/C 240		
2638A/60/C 220		
2638A/60/C 240		
2638A/20 220-6/CN		
2638A/40 220-6/CN		
2638A/20 220		
2638A/20 240		
2638A/40 220		
2638A/40 240		
2638A/60 220		
2638A/60 240		
1586A 240/C		
1586A/1DS 240/C		
1586A/2DS 240/C		
1586A/1HC 240/C		
1586A/2HC 240/C		
1586A/DSHC 240/C		